

## VT 1200S Customer Training

Introduction & Sectioning Technique

2025/11/06

龐德生技 Henry



# Agenda

Advancing Cancer Diagnostics  
Improving Lives



VT1200S Introduction

Hardware, operation & maintenance

Sectioning technique

Hands-on practice

# VT1200S Introduction

# What is Vibratome?

Advancing Cancer Diagnostics  
Improving Lives



- A Vibratome is a microtome with a vibrating blade to section soft tissue under physiological conditions
- Special sectioning technique: the blade vibrates laterally and moves slowly forward through the specimen to reduction of cutting forces
- Limitation in section thickness compared to cryostats or rotary microtomes
  - Fixed tissue – down to 10µm
  - Unfixed
  - Fresh tissue mainly 100 – 300µm

# Why using a Vibratome?

Advancing Cancer Diagnostics  
Improving Lives



Advantage of sectioning under physiological conditions

- No thermal effect
- No chemical treatment necessary
- No embedding necessary
- Chilled physiological buffer surrounds the tissue while sectioning: Antigen protection / reduction of enzyme activity

# VT1200S – Intended use & Features

Advancing Cancer Diagnostics  
Improving Lives



- RESEACH USE ONLY! Not for diagnosis!
- For SOFT fixed / unfixed tissue
- Knife holder technology: magnetic
- Adjustable amplitude up to 3mm
- Fixed frequency (85 Hz)
- Clearance angle adjustable: 15 °, 18 ° and 21 °
- Sectioning Mode: semiautomated and automated (single and continuous)
- Vibrocheck technology




# VT1200S – Application

Advancing Cancer Diagnostics  
Improving Lives



## Soft fixed / unfixed tissue in buffer

- Brain tissue
  - Fixed or unfixed
  - Electrophysiology / patch clamping
    - Rather thick (min. 30µm) sections needed
    - Vibrocheck technology ensures that the blades moves absolutely parallel through the tissue
    - Living cells, even on section surfaces
  - Fluorescence work
  - Immuno-staining
  - Cell cultures
- Other tissue
  - Mainly cell cultures
    - Liver, Heart, Kidney



# VT1200S – Application

Advancing Cancer Diagnostics  
Improving Lives



- Vibratomes are not designed to produce “thin” sections
  - The advantage of avoiding embedding / freezing however limits the “thinness” of the section – for living tissue to record cells – too thin sections would result in dead cells (cell size approx. 30µm)
  - If you are looking for thin sections below 10µm – we recommend:
    - rotary microtomes (embedded tissue)
    - cryostats (frozen tissue)
    - sliding microtomes with freezing attachment (fixed – frozen tissue)



# VT 1200S Hardware & operation

# VT1200S – Hardware overview


Advancing Cancer Diagnostics  
Improving Lives



Socket (right hand side)



Cutting head



Blade holder

Control panel

Buffer tray and ice bath

Magnifier (add-on)




Microscopy (add-on)

Adjustable 2-arm LED light (add-on)



Power socket & power switch




# Hardware – Cutting head & blade holder


Advancing Cancer Diagnostics  
Improving Lives




- Magnetic knife holder system
- Specimen retraction
- Usable for razor: injector /sapphire blades
- Clearance angle adjustable: 15 ° , 18 ° and 21 °
  - Recommended setting: 18 °
  - For 15°, the effective clearance angle is "0"



Hole for tile up  
blade holder



15 ° 18 ° 21 °



Hole for remove  
the blade holder



# Hardware – Insert a blade


Advancing Cancer Diagnostics  
Improving Lives

**Leica**  
BIO SYSTEMS

1. rotate it 90°  
clockwise (tile up) with  
allen key




2. insert allen key from  
the top through the  
opening (10) into the  
blade holder (11) and  
open the blade holder



3. Hold the entire razor  
blade (not separated) on  
the left and right with  
both hands and insert it  
into the blade holder.  
Hook in the blade over  
bottom pressure plate



4. insert allen key from  
the top through the  
opening (10) into the  
blade holder (11) and  
tighten the blade  
holder




# Hardware – Specimen plate, Buffer tray & Ice tray

Metal buffer tray  
(cannot be autoclaved)

Plastic buffer tray  
(can be autoclaved)

Lid



Ice tray



Manipulator and  
specimen plate



# Hardware – Install ice tray & buffer tray

Advancing Cancer Diagnostics  
Improving Lives



1. Insert buffer tray onto ice tray
2. Cover the buffer tray with the acrylic lid
3. fill the ice tray with crushed ice
4. Uncover the buffer tray and fill it with cooled buffer solution
5. Pulled the ice tray forward for clamping on device
6. Clamp it down by pushing the lever towards the rear

# Hardware – Control panel

Advancing Cancer Diagnostics  
Improving Lives



In semiautomatic mode, the desired section thickness is selected with Knob2. Values appears in  $\mu\text{m}$ . After each sectioning process has been completed, the  $\mu\text{m}$  display is set to "0"

**SPEED:**  
blade feed rate  
0.01 to 1.5 mm/s

**AMPLITUDE:**  
0 to 3 mm

**\*AUTO FEED:**  
max. 1000  $\mu\text{m}$

**LIGHT  
ON/OFF**

**ROTARY KNOB 1:**  
For speed,  
amplitude, auto feed

**ROTARY KNOB 2:**  
For moving  
specimen receptacle  
vertically

**STEP SIZE (use  
with knob 2)**  
Choose 1, 10 or  
100  $\mu\text{m}$



**CLEAR:**  
Display for the  
current specimen  
holder position /  
reset the display  
to 0 at any point

**MENU:**  
User profile

# Hardware – Control panel

Advancing Cancer Diagnostics  
Improving Lives



## UP:

Press and hold: specimen receptacle moves upward  
Release: stop moving  
lights up when pos. is the highest

## AUTO / MAN:

Switch between auto and manual (semiauto) sectioning

## SINGLE / CONT:

Switch between single and continuous sectioning

## PAUSE:

Stop immediately  
Press again to restart the current process  
Press RUN / STOP to discontinues the sectioning process

## \*CUTTING WINDOW:

Set the cutting path edge  
Holding the key down for 3s to activate / deactivated  
LED lights up when edge is set

## DOWN:


Press x1: specimen receptacle auto moves to lowest pos.  
lights up when pos. is the lowest  
Press again: stop moving

## BLADE FORWARD / BACK:

Press and hold: blade moves forward / back  
LED lights up when reach the end point

## RUN / STOP:

Press to start  
In CONT mode, blade will stop after the current sectioning is done





# Operation – Retraction between auto & semiauto sectioning

Advancing Cancer Diagnostics  
Improving Lives



## What is retraction?

- specimen is lowered to the desired retraction value after each cut to prevent the specimen surface and the blade from coming into contact while the blade is being retracted.
- **Recommended for section soft tissue**
- **In semiautomatic sectioning mode**, a manual feed to the desired section thickness must be carried out before each cut. There is no automatic specimen retraction; however, retraction can be performed manually.
- **In automatic mode**, an automatic feed (AUTO FEED) to the selected section thickness is carried out before each cut, and the specimen is lowered to the desired retraction value after each cut.

# Operation – Automatic sectioning

Advancing Cancer Diagnostics  
Improving Lives



1. Specimen preparation (Put specimen in buffer tray in ice tray)
2. Switch on VT1200S
3. Tile up blade holder to insert the blade,
4. Load the ice tray (with specimen inside) to VT
5. Tile down blade holder and to the desired clearance angle
6. Choose auto mode & single / cont, set speed, amplitude, auto feed thickness
7. According to the specimen position, set the cutting window and specimen receptacle position
8. Start the run
9. When finish sectioning, lower the specimen receptacle, replace with another specimen plate with the manipulator or remove the whole ice tray
10. Remove the blade



# Operation – Semiautomatic sectioning

Advancing Cancer Diagnostics  
Improving Lives



1. Specimen preparation (Put specimen in buffer tray in ice tray)
2. Switch on VT1200S
3. Tile up blade holder to insert the blade,
4. Load the ice tray (with specimen inside) to VT
5. Tile down blade holder and to the desired clearance angle
6. Choose semiauto mode, set speed, amplitude
7. According to the specimen position, set specimen receptacle position (=thickness) and blade position
8. Start the run
9. Lower a bit the specimen receptacle with knob 2, move the blade back, move up the specimen receptacle ( $\mu\text{m}$  display is set to "0") (manual retraction)
10. Repeat the step 7 to 10
11. When finish sectioning, lower the specimen receptacle, replace with another specimen plate with the manipulator or remove the whole ice tray
12. Remove the blade

# Maintenance

Advancing Cancer Diagnostics  
Improving Lives



Before each cleaning, carry out the following preparatory steps:

1. Switch off the main switch on the side of the instrument.
2. Place the magnifier cover on the magnifier.
3. Remove the blade from the blade holder and dispose of it safely.
4. Pull the ice tray and buffer tray off of the dovetail guide and place them on the stage.
5. Remove and empty out the buffer tray. Dispose of the contents of the buffer tray properly.
6. Remove the specimen plate and lay it flat on the stage.
7. Remove the specimen using a single-sided blade and remove any cyanoacrylate adhesive residue from the specimen plate.

# Maintenance

Advancing Cancer Diagnostics  
Improving Lives



- Instrument & outside surfaces
  - Cleaned with a mild commercial household cleaner or soap water (DO NOT USE SOLVENT CONTAINING XYLENE / ACETONE) and then be wiped with a cloth.
  - completely dry before use
- Specimen plate / buffer tray / ice tray / Allen key / manipulator
  - Wash thoroughly with de-ionized water and finally with 100% Ethanol
  - Dry the parts on a tissue paper before re-installing them
- Cleaning the blade holder
  - Remove and spray with alcohol. Wiped down with a piece of cellulose and placed on a cellulose towel to dry completely
- Cleaning the sapphire blades
  - using an alcohol-based solution or acetone



# Sectioning technique

# Sectioning technique

Advancing Cancer Diagnostics  
Improving Lives



- Sectioning speed
  - select a moderate knife travel forward speed to avoid that the specimen is pushed while sectioning
  - The “Rule of thumb” is that softer tissues require slower forward speed, higher amplitude and frequency, steeper blade angle and cold buffer
  - Fresh brain tissue: amplitude 1.5 - 1.8 mm / speed: very slow
  - Fixed brain tissue: amplitude 0.8 –1.0 mm / speed: slow
  - Liver: amplitude 3mm / speed: 10µm/s –the slowest speed needs to be used
- Sectioning under physiological conditions
  - work with cold buffer
    - Firm specimen is easier to be cut
    - Crushed ice can be used to keep the buffer cooled or more elegantly, a double walled buffer with a circulation Julabo Chiller device
- When sectioning fresh tissue,
  - oxygenate the tissue while sectioning to ensure viability of the cells



# Sectioning technique

Advancing Cancer Diagnostics  
Improving Lives



- Specimen preparation: Tape-mounted sectioning




# Sectioning technique

Advancing Cancer Diagnostics  
Improving Lives



- Specimen preparation: Agarose embedded block
  - For very small and soft tissue



Recommended agarose with low melting point of 37 to avoid heating artifact

*For Research Use Only. Not For Use In Diagnostic Procedures.*


Leica Biosystems Proprietary Information

# Sectioning technique

Advancing Cancer Diagnostics  
Improving Lives

**Leica**  
BIOSYSTEMS

- Specimen preparation: Agarose barrier





# Thank you

